

Fixing and Maintenance Instructions

Commodes

N.B. This equipment must be installed by a competent installer/ technician.

Instructions for:

Kent Commode	154kg (25st.) weight limit
Essex Commode	154kg (25st.) weight limit
Sussex Commode	154kg (25st.) weight limit
Norfolk Commode	154kg (25st.) weight limit
Solo Skandia	127kg (20st.) weight limit
President Toilet Seat and Frame	up to 154 kg (25st) weight limit
Ashby and Cosby Lux	154kg (25st.) weight limit
Grand Cosby Lux	248kg (39st.) weight limit

Do not exceed the weight limits stated – doing so could put the user at **risk**.

1. All Commodes

Please examine your commode for any visible damage before use. If you see any damage or suspect a fault, please contact your supplier.

In this case, please do not use your product as it may compromise your safety.

2. Adjustment of legs (if applicable)

Remove the 'E' clip from the extension leg and extend all legs evenly to the desired height ensuring you do not exceed the last punched hole of the extension leg. (see fig. 1)

Once you have achieved the height you need, re-install the 'E' clip ensuring that it has passed through the extension leg and is fitting neatly and snugly around the commode leg. (see fig. 2)

Make sure that all legs are sitting evenly onto the floor; **never** extend legs in a manner that the commode is slanting. **This will compromise the safety of the user.**

3. Surrey Commodes only

The Surrey Commode has detachable arms; to remove these simply lift the arms out the locating tubes welded to the side of the commode frame. Always ensure that the arms are re-fitted neatly and evenly into the locating tubes before use. Failure to do so will compromise the safety of the user.

4. Commode Buckets (Kent, Essex, Sussex, Norfolk and Solo Skandia Commodes)

To remove the commode bucket, simply lift the toilet seat and fix the lid provided onto the bucket to avoid splashing and carefully raise the bucket and dispose of the contents into a proper sluice or WC. Never use an open drain, sink or basin. This could cause a serious health risk. Clean out the residue with a suitable disinfectant and refit to the commode after use

5. Cleaning (All Commodes)

All commodes must be cleaned using a non-abrasive cleaner or can be pressure washed. However, always ensure that any chemicals you use with your pressure washing will not affect chrome plate, polymer coating or epoxy powder coating. After cleaning use a water dispersing agent such as WD40 on all joints and castors to reduce the chance of premature rust or seizure.

6. Reissue

If you reissue or are about to reissue this product, please thoroughly check all components for their safety. This includes;

P.A.T. test for electrical items, tightness of all nuts/ bolts/ castors and other screw in/ bolt in/ push in components. Also check all upholstery for security, splits, etc.

If in any doubt, please do not issue or use, but immediately contact your supplier for service support.

